

ELECTROLYTES REPLENISHMENT: Overview

The Equicure Electrolytes Replenishment is an injector providing fast recuperation of the horse during and after exertion.

Its composition gives fast mineral salts replenishment to counterbalance the mineral deficiencies caused by exertion and/or to give back a "fillip" to the horse.

Its formulation is based on products made from the transformation of plants selected for their electrolytic contribution, such as valerian, chamomile or yellow jasmine.



What are the advantages of the Equicure Electrolytes Replenishment ?

- *An easy to use product:* the oral syringe packaging means you can carry it with you during hacks, competition,... It is immediately usable, designed to be easily swallowed by the animal, without dilution, nor particular manipulation.
- *A product that can be used for all:* use during gestation, foaling, nursing and for very young horses.
- *A custom-made product:* the quantity of electrolytes by injector has been optimized to avoid absorption losses and to maximize the performance of your horse.
- *A non-doping product:* the Equicure Electrolytes Replenishment doesn't show up positive in the animal during a doping test and can therefore be used in competition.

When should the Equicure Electrolytes Replenishment be used?

Electrolytes Replenishment is used in the following situations:

- During exercise, competition, work or all other activity causing sweating and requiring an instantaneous electrolytes replenishment and fast energy renewal - enabling the horse to continue the activity in good condition.
- After exertion resulting in sweating (work, race, transportation, ...) requiring a fast recuperation in electrolytes.

What are the properties of the Equicure Electrolytes Replenishment?

The Equicure Electrolyte Replenishment allows:


- A faster recuperation of your horse during and / or after exertion.
- Stay fit longer during exertion.
- To give back a "fillip" to your horse and will perform much longer.
- To be able to alternate between training / competition more easily.
- To quickly restore the loss of electrolytes caused by profuse sweating.

How should the Equicure Electrolyte Replenishment be used?

Place the injector at the corner of the animal's mouth and press the piston to empty its content.

The recommended doses are:

- Horses : 1 / 2 to 1 injector
- Ponies – foals: 1/ 2 an injector



Distributed by :